

**ST. JOHN HOSPITAL AND MEDICAL CENTER
LAB COAT ORDER FORM**

Three (3) lab coats are provided to SJH&MC residents at the start of their program. Coats are stored in the Linen Room.

MEASURING

To get an accurate chest measurement, place measuring tape around fullest part of chest, keeping it up under arms and across the shoulder blades. The measurement in inches is your size. **Be sure measurement is made over the clothing to be worn. Please be as accurate as possible, coats cannot be returned for corrections to sizing.**

Straight back, side pockets are available in regular or long lengths (long adds 2" to sleeves and 4" overall length).

Sizes (Regular or Long lengths):

- | | |
|--------------------------|--------------------------|
| <input type="radio"/> 30 | <input type="radio"/> 44 |
| <input type="radio"/> 32 | <input type="radio"/> 46 |
| <input type="radio"/> 34 | <input type="radio"/> 48 |
| <input type="radio"/> 36 | <input type="radio"/> 50 |
| <input type="radio"/> 38 | <input type="radio"/> 52 |
| <input type="radio"/> 40 | <input type="radio"/> 54 |
| <input type="radio"/> 42 | |

ORDER:

- **NAME (please print)** _____
- **LENGTH (circle one)** REGULAR or LONG
- **QUANTITY:** 3 _____

RETURN THIS FORM WITH ALL OTHERS TO THE GME OFFICE